

La gestione delle prove dirette per Brucella nei piani di eradicazione

Antonio Battisti

Direzione Operativa Diagnostica Generale
IZS Lazio e Toscana, Sede Centrale di Roma

Teramo, 27 febbraio 2015

Direzione Operativa Diagnostica Generale

- Struttura Complessa presso la Sede Centrale di Roma
- Centro di riferimento (bi)-regionale per agenti zoonosici Speciali (inclusi alcuni agenti di classe di rischio 3: Brucella spp., M. bovis)
- Gestisce prove dirette per Brucella per l'IZSLT

Capacità di Laboratorio

- Riferimenti per le prove colturali e molecolari:
**OIE Terrestrial Manual 2012. BOVINE BRUCELLOSIS.
CH. 2.4.3. (1.b; 1.c; 1.d)**

Version adopted by the World Assembly of Delegates of the OIE in May 2009.

- La Struttura esegue prove colturali dirette per Brucella secondo POS accreditata (ISO/IEC 17025 Standard)
- Validazione interna (Manual OIE, Ch 1.1.15 OIE);
- In seguito, evidenza di mantenimento validità con circuito interlab CRN, IZS Abruzzo e Molise

Approccio alla validazione del metodo

- Validazione secondo principi ISO & OIE Manual, fin dal 2003.

Attualmente:

- Ch. 1. 1.15: Principles and methods of validation of diagnostic assays for infectious diseases (2013)

Validazione del metodo

- **Valutazione caratteristiche performance della prova (D-Se & D-Sp):** Preparazione di campioni POSITIVI e NEGATIVI di Riferimento (n=30 e n=90) nel 2003
- Tabella di contingenza 2X2
- D-Se, D-Sp (e 95% I. C.)
- Accettabilità: D-Se: $\geq 95\%$; D-Sp $\geq 98\%$

1.b Cultures


- Basal media: Brucella medium base+5% equine serum
- Selective media:
 - Farrell's Medium (VCN supplement)

Farrell's medium (Farrell, 1974), which is prepared by the addition of six antibiotics to a basal medium. The following quantities are added to 1 litre of agar: polymyxin B sulphate (5000 units = 5 mg); bacitracin (25,000 units = 25 mg); natamycin (50 mg); **nalidixic acid** (5 mg); nystatin (100,000 units); vancomycin (20 mg).

Incubation: 8-10 days

(Feti e campioni da aborti: alte concentrazioni di Brucella spp.)

Colture primarie (a 7 gg) di feto ovino abortito, 12/2014,
su Farrell's medium (*B. melitensis* 3)



Alternative selective solid medium

-Modified Thayer-Martin's medium (for some *B. abortus* & *B. melitensis* biotypes)

GC medium base, Hb 10g/L, colistin (7.5 mg/litre), vancomycin (3 mg/litre), nitrofurantoin (10 mg/litre), nystatin (100,000 I U/ litre = 17.7 mg) and amphotericin B (2.5 mg/litre)

Selective enrichment

-Brucella broth/Tryptose (Trypticase)-soy broth+5% equine serum for enrichment:

supplemented with an antibiotic mixture of **at least amphotericin B (1 µg/ml), and vancomycin (20 µg/ml)** (all final concentrations).

The enrichment medium should be incubated at 37°C in air supplemented with 5–10% (v/v) CO₂ **for up to 6 weeks**, with **weekly subcultures** onto solid selective medium

Cultures (OIE Manual of Standards, 2009)

- Singoli reattori, casi sospetti o dubbi:
- Media di arricchimento incubati at 37°C in 5-10% (v/v) CO₂ fino a 6 settimane (operativamente di solito 3), con:


In parallelo

-subcolture media selettivi solidi

-PCR testing

Entrambe ad intervalli settimanali

Caratteristiche fenotipiche macroscopiche e microscopiche di *Brucella* spp.



Identificazione degli isolati a livello di Genere

1.c Identification and typing

- Test biochimici differenziali (min. CAT, OX, URE), oppure Gallerie di test biochimici differenziali (es. API 20NE)
- slide agglutination: acriflavina, siero polivalente, o anche monovalenti A e M)
- **PCR testing** (Brucella spp.):
 - **End-point PCR:**
 - specifici primers per amplificazione di sequenza di 223 bp di antigene immunogenico di 31 kDa di Brucella abortus (**BCPS31**, Baily et al., 1992; primers: Elfaki et al., 2005)
 - **Real-Time PCR** (Bounaadja et al., 2009)


Nucleic acid recognition methods...

- **1. d) Nucleic acid recognition methods**

“The **PCR**, including the **real-time** format, provides an additional means of detection and identification of *Brucella* sp. (**Bricker, 2002; Bricker et al., 2003; Bricker & Halling, 1994; 1995; Garcia-Yoldi et al., 2006; Hinić et al., 2008; Ocampo-Sosa et al., 2005**). (....)”

Capitolo Brucellosi dell'OIE Manual non aggiornato dal 2009:
non cita la Ref. Real-Time PCR del EU-RL/OIE-RL

BCPS31 PCR (223 bp amplicon)



Da matrici con bassa concentrazione di *Brucella* spp.

Altri approcci molecolari

- Per l'evidenziazione di sequenze nucleotidiche specifiche per *Brucella spp.* da colture (arricchimento) e materiali biologici (es. Latte, tessuti mammari, linfatici):

Real-Time PCR basata su **IS 711 locus**

Detection & ID by Real-Time PCR

Veterinary Microbiology 137 (2009) 156–164



Contents lists available at ScienceDirect

Veterinary Microbiology

journal homepage: www.elsevier.com/locate/vetmic



Real-time PCR for identification of *Brucella* spp.: A comparative study of IS711, *bcs*31 and *per* target genes

Lotfi Bounaadja^{a,b}, David Albert^c, Benoît Chénaïs^a, Sylvie Hénault^c,
Michel S. Zygmunt^d, Sylvie Poliak^b, Bruno Garin-Bastuji^{c,*}

^a Université du Maine, Laboratoire de Biologie et Génétique Evolutive, EA2160 Mer Molécules Santé, 72085 Le Mans, France

^b Laboratoire Départemental de la Sarthe, 72000 Le Mans, France

^c OIE/FAO and EU Community Reference Laboratory for Brucellosis, French Food Safety Agency (AFSSA), 23 avenue du Général de Gaulle, 94706 Maisons-Alfort Cedex, France

^d UR1282, Infectiologie Animale et Santé Publique (IASP), Institut National de la Recherche Agronomique (INRA), 37380 Nouzilly, France

ARTICLE INFO

Article history:

Received 18 July 2008

Received in revised form 22 December 2008

Accepted 29 December 2008

Keywords:

Brucella

Real-time PCR

*bcs*31

per

IS711

Genus-specific identification

ABSTRACT

Culture is considered as the reference standard assay for diagnosis of *Brucella* spp. in humans and animals but it is time-consuming and hazardous. In this study, we evaluated the performances of newly designed real-time PCR assays using TaqMan[®] probes and targeting the 3 following specific genes: (i) the insertion sequence IS711, (ii) *bcs*31 and (iii) *per* genes for the detection of *Brucella* at genus level. The real-time PCR assays were compared to previously described conventional PCR assays targeting the same genes. The genus-specificity was evaluated on 26 *Brucella* strains, including all species and biovars. The analytical specificity was evaluated on a collection of 68 clinically relevant, phylogenetically related or serologically cross-reacting micro-organisms. The analytical sensitivity was assessed using decreasing DNA quantities of *Brucella ovis*, *B. melitensis* bv. 1, *B. abortus* bv. 1 and *B. canis* reference strains. Finally, intra-assay repeatability and inter-assay reproducibility were assessed. All *Brucella* species DNA were amplified in the three tests. However, the earliest signal was observed with the IS711 real-time PCR, where it varied according to the IS711 copy number. No cross-reactivity was observed in all three tests. Real-time PCR was always more sensitive than conventional PCR assays. The real-time PCR assay targeting IS711 presented an identical or a greater sensitivity than the two other tests. In all cases, the variability was very low. In conclusion, real-time PCR assays are easy-to-use, produce results faster than conventional PCR systems while reducing DNA contamination risks. The IS711-based real-time PCR assay is specific and highly sensitive and appears as an efficient and reproducible method for the rapid and safe detection of the genus *Brucella*.

© 2008 Elsevier B.V. All rights reserved.

Parametri per la Validazione (Real-Time) PCR


- I parametri considerati nel processo di validazione del metodo di PCR Real Time per la ricerca del gene target IS711 in *Brucella* spp. sono:
- Specificità analitica (Analytical Specificity, **Asp**, Ch. 3.6.3, 1.2)
- Limite di rivelazione (LOD) o Sensibilità analitica (Analytical Sensitivity, **Ase**, Ch. 3.6.3, 1.3)
- Robustezza
- Sensibilità Diagnostica
- Specificità Diagnostica

Tipizzazione a livello di specie e biotipo

- Attualmente la DO Diagnostica Generale non investe in risorse per la successiva caratterizzazione.
- Isolati inviati a CRN Brucellosi, IZS Abruzzo e Molise, Teramo

Brucella isolates, Latium

BRUCELLA ISOLATES LATIUM 2
by BRUCELLA SPECIES BIOTYPE



Brucella suis

- 5 biovars, di cui solo le prime tre ritenute infettanti il suino
- In Europa *B. suis* 2 è descritta in alcune popolazioni di cinghiali dell'Europa dell'Est e anche in MS (Italia, Portogallo, Francia)
- La situazione epidemiologica non è chiara nella maggior parte dei MS, ma si sospetta che molte popolazioni di cinghiali mantengano l'agente in EU



Allevamento suino semi-brado, *B. suis* 2









Contents lists available at SciVerse ScienceDirect

Veterinary Microbiology

journal homepage: www.elsevier.com/locate/vetmic



Short communication

The presence of *Brucella ceti* ST26 in a striped dolphin (*Stenella coeruleoalba*) with meningoencephalitis from the Mediterranean Sea

Patricia Alba^a, Giuliana Terracciano^a, Alessia Franco^a, Serena Lorenzetti^a, Cristiano Cocumelli^a, Gianluca Fichi^a, Claudia Eleni^a, Michel S. Zygmunt^{b,c}, Axel Cloeckaert^{b,c}, Antonio Battisti^{a,*}

^a Istituto Zooprofilattico Sperimentale delle Regioni Lazio e Toscana, Via Appia Nuova 1411, 00178, Rome, Italy

^b INRA, UMR1282 Infectiologie et Santé Publique, F-37380 Nouzilly, France

^c Université François Rabelais de Tours, UMR1282 Infectiologie et Santé Publique, F-37000 Tours, France



ARTICLE INFO

Article history:

Received 7 November 2012

Received in revised form 18 January 2013

Accepted 22 January 2013

Keywords:

Brucella
Brucella ceti dolphin type
 Cetacea
 Striped dolphin
 Infectious diseases
 Marine mammals
 Mediterranean Sea
 Molecular characterization
 Pathology
Salmonella
Toxoplasma gondii
 Zoonoses

ABSTRACT

Brucella spp. was isolated from brain, lung and intestinal lymph nodes of a dead adult male striped dolphin (*Stenella coeruleoalba*) found stranded on the Tyrrhenian coast (Tuscany, Italy) of the Mediterranean Sea in February 2012. *Brucella* spp. was associated with moderate to severe lesions of meningoencephalitis. A co-infection by *Toxoplasma gondii* was also demonstrated at brain level by means of molecular and histopathologic methods. The *Brucella* isolate was further characterized based on a fragment-specific polymerase chain reaction (PCR) approach, consisting of a set of five specific PCRs, targeting specific chromosomal IS711 locations for marine mammal *Brucellae*, as described previously. The isolate was thus classified as *Brucella ceti* 1; V fragment-positive (or *B. ceti* dolphin type), according to previous studies. Multi Locus Sequence Analysis demonstrated that the isolate belongs to Sequence Type 26, while *omp2* (*omp2a* and *omp2b* genes) sequence analysis further confirmed the isolate belonged to this group of strains. This is the first report of *Brucella* spp. from marine mammals in the Mediterranean Sea, and represents a further observation that this strain group is associated with hosts of the Family *Delphinidae*, and particularly with the striped dolphins, also in the Mediterranean area, thus constituting a further biological hazard of concern for this vulnerable subpopulation. © 2013 Elsevier B.V. All rights reserved.

Grazie dell'attenzione

Un grazie particolare al personale delle D. O. Diagnostica Generale

